

MW40



# Obstetric Training Simulator -complete set-

**Comprehensive maternal body set for skills training  
in vaginal examinations during labor and delivery**



## FEATURES

- The external genitalia is made of innovative material that excels in elasticity, durability and resilience for effective and repetitive training.

## SKILLS

- Cervical examination
- Delivery assistance
- Perineorrhaphy (perineal repair)


## KEY FEATURES

The external genitalia is made of innovative material that excels in elasticity, durability and resilience for effective and repetitive training.



## LANDMARKS

- velamen (amniotic sac)
- vaginal canal
- ischial spine
- pubic symphysis
- coccyx
- anterior fontanelle



**KEY FEATURES**

**Cervical examination**

5 types of cervical dilation models and fetal descent positions can be interchanged with a single action.

- | Cervical scoring in the first stage of labor
  - Determination of station using landmarks of the ischial spine and the pubic joint to assess the descent of the body
  - Assessment of cervical dilation
- | Speculum examination
- | Observation and cervical cytology
- | Urinary catheterization



**Delivery assistance**

**Smooth and realistic delivery**



- | Delivery assistance for different delivery positions (dorsal, all-fours)
- | Perineal protection
- | Delivery of fetus (forceps delivery, vacuum extraction)
- | Clamping, tying and cutting of the umbilical cord
- | Delivery of the placenta
- | Insertion of urinary catheter



**Perineorrhaphy (perineal repair)**

Repetitive suture training for the first-degree perineal tears.

**First-degree perineal tear repair**

- | Single interrupted suture
- | Vertical mattress suture



**DESCRIPTIONS**

**SPECIFICATIONS**

Size: W74 x D44 x H33.5cm / W29.1x D17.3 x H13.2 inch Weight: 11.2kg / 24.7lbs

**MATERIALS**

Soft resin, Latex free

**SET INCLUDES**

- 1 maternal body
- 1 genitalia unit (for pelvic examination and delivery assistance)
- 1 talcum powder
- 1 lubricant (500ml)
- 1 towel
- 1 instruction manual

**Parts for pelvic examination:**

- 1 vaginal unit
- 1 set of cervical dilation inserts (5 variations)
- 1 holder for the cervical dilation inserts

**Parts for Delivery Assistance:**

- 1 fetus
- 1 placenta
- 3 umbilical cords for omphalotomy
- 1 tray
- 5 velamen sheets

**Parts for Perineorrhaphy:**

- 1 genitalia unit (for perineorrhaphy)